


Before you can launch your third-party imaging software, you must first set it up to work within the ABELDent program. Keep in mind the following when setting up:

- Set up the imaging **integration** on all computers to be used for image capturing.
- See **Third Party Imaging Packages** for a complete list of the bridges ABELDent has developed.

1. To Set Up Third Party Imaging Software:

1. On the **Imaging** tab right-click and choose **Setup Integration**.

The **Integration Selection** dialog appears.



2. Click on the product you want to bridge to and click **Enable**.

A window will appear and prompt you for the product's .exe file.

3. Locate and select the **.exe** file and click **Open**.

4. Follow the on-screen prompts.

5. Once the **integration** setup has been completed you can launch your imaging software from either the patient's imaging tab, or by clicking the **Launch Integration** icon on the toolbar.