

ABELDent Advantage Cloud Hardware Requirements

Important Note: If you are purchasing new hardware, we recommend that you purchase hardware that meets or exceeds our recommended specifications.

****Don't rely entirely on paper copies of these specifications. Always check our web site for the most current requirements.**

| Component | Recommended | Details |
|------------------------------------|--|--|
| Server is not required for cloud | | A local server is not required for ABELDent Advantage running in the cloud. Your radiography or imaging solution may require a server. Please check with your imaging vendor for details on their requirements. |
| Client workstations | Intel i5 2.5 GHz or better, 4GB RAM, 120GB Hard disk(SSD), DVD Windows 10 Pro We recommend processors rated 5000 or above on CPU Benchmarks https://www.cpubenchmark.net/high_end_cpus.html | ABELDent Advantage is designed for use with Windows 10 Professional (English only). Most current radiographic and imaging programs also work with Windows 10, please confirm support for your version with your imaging vendor. ABELDent Advantage does not support older or home operating systems. Use of unsupported operating systems can lead to failure of the software to operate or damage to system data. Recent Dual or Quad core processors will perform well. These include Core processors, i5, i7 from Intel and AMD equivalents. |
| Display | 21" with 1920x1080 resolution | |
| Backup Unit | Not Required for ABELDent Advantage data | ABELDent Advantage Cloud is a hosted offering so we take care of backups. Your other radiography, imaging or other applications may have data that you need to protect so you may require an Internet backup service or external drives (encrypted) to protect that data. All external media, whether for backup or other purpose, should be encrypted before copying confidential information to it. |
| Network Interface Card | 1 Gbps NIC | In all situations where computers share information on a network, each PC will require a network interface card. |
| Switch | 1Gbps switch | In a typical small business network the cables converge at an Ethernet switch. Gigabit speed switches are recommended. A port is required for each PC. Note: Wired connections are preferred as some wireless network connections are not reliable enough for database connections. |
| Cable | Cat 6a, 6 or 5e Ethernet (See Details column) | Ethernet cable will be required in your office running from each computer to a switch. The switch is normally located centrally in your office, possibly in a wiring closet. |
| Internet Connectivity | High speed Internet is required | Low to midrange Business plans are suitable for small practices. Practices with 6 or more workstations should consider higher plans. |
| Laser Printer | Any printer with appropriate Windows Drivers (32MB RAM recommended) | All recent revisions of ABELDent support Windows GDI printing. HP P1606DN, 2000, 2015, 3005, and 4000 series printers are known to work with ABELDent. |
| Optional Color Inkjet Printer | Any color inkjet with appropriate Windows Drivers | For printing Charts & Images. Check the number of pages per minute the printer will print. Some are quite slow. |
| Uninterruptible Power Supply (UPS) | Optional as the server in the cloud is protected from power outages. | If you want to maintain connectivity to the cloud during local power interruptions, you can protect your modem/router switching equipment and workstations of your choice with UPS(s). |
| Imaging | USB Camera | Most modern cameras plug into a USB port and do not require Video capture cards. |
| Optional Scanner | TWAIN | Any TWAIN compliant scanner will work with ABELDent |
| Optional Signature Pads | ePad (VP9801) or ePad-ink (VP9805) | ABELDent Advantage only support these two models. They can be ordered online from the ePadLink website at http://www.epadlink.com/ |
| Sound Card & Speakers | See Details column. | A sound card and speakers or headphones are required for the ABELDent training videos, Optio Patient Education, and to view demo videos. |
| Headset | See Details column. | A good quality USB headset will be required for users who will be using voice activated Perio Charting. |
| Antivirus Software | See Details column. | AntiVirus/AntiMalware/AntiRansomware Software, configured to automatically check for definition updates daily, is strongly recommended for every PC connected to the Internet. |
| Firewall / Router | See Details column. | A business grade router/security appliance is strongly recommended for all Networks with a high-speed Internet connection. Note: Wired connections are preferred as some wireless network connections are not reliable enough for database connections |